# Ebulliometer

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**Notice d'utilisation  
Instruction manual  
Manual de uso  
Version 2012  
Ref. 160350T/160350D**

**LABORATOIRES DUJARDIN-SALLERON**

## LABORATOIRES DUJARDIN-SALLERON

872, route de la gare  
37210 NOIZAY - FRANCE

Tél. : +33 (0)2.47.25.58.25  
Fax : +33 (0)2.47.25.58.30

[info@dujardin-salleron.com](mailto:info@dujardin-salleron.com)  
[www.dujardin-salleron.com](http://www.dujardin-salleron.com)

# SOMMAIRE

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

# 1- General description of the instrument and the accessories

## Instrument :

1. Introduction of samples
2. Condenser upper outlet
3. Condenser Réf. 160350/1
4. Cold water inlet
5. According to the model:
  - Mercury thermometer Réf. 690117E
  - Temperature probe Réf. 160350/2
6. Thermometer stopper Réf. 160350/3
7. Boiler Réf. 160350/1
8. Heating cartridge Réf. 160350/4
9. Drain valve
10. LCD screen for temperature  
(Model Ref.160350D)
11. Heating switch
12. 6-plug DIN connector – connection of the  
temperature probe (Model Ref.160350D)
13. Adapter connector - Power supply 24V DCa
14. ON/OFF switch

## Accessories:

- |                                                     |               |
|-----------------------------------------------------|---------------|
| 15. Disc for electric ebulliometer                  | Réf. 160350/5 |
| 16. Standard wine for calibration                   | Réf. 100550   |
| 17. Defoamer for ebulliometer                       | Réf. 990001   |
| 18. Chromosulfuric acid<br>(Supplied conditionally) | Réf. 972500   |
| 19. Power cable                                     |               |
| 20. Cooling system and drain hoses                  |               |
| 21. Dropping-pipette with bulb                      | Réf. 630011   |

## 2 - Setting up of the instrument

### • Installation of the instrument

- Put the instrument on a flat surface
- Tightly and carefully push in the boiler on the cartridge holder. You may add some water on the washers to make this operation easier.
- Connect the instrument on the power supply.
- Brancher l'appareil sur secteur.

*Installation of boiler on stand*


### • Installation of the probe

Place the probe on the boiler without modifying the position of the silicon gasket on the probe, adding some water on the silicon gasket if necessary in order to push it in at the maximum. The end of the probe must be situated between 0,5 and 1 cm above the end of the heating element.

In case of the model with mercury thermometer, place it the same way, making sure that the mercury bulb or the probe is 1 cm above the heating cartridge.

In both case, the measuring stem must be on the same axis than the heating cartridge.

### • Installation of the cooling system and tubing

Attach the water connections to the condenser.

- Cold water inlet to the bottom;
- Water outlet to the top.

Attach additional tubing to the boiler drain and place the tubing into a sink to drain easier.

## 3 - Calibration

The calibration allows to measure the boiling point of a standard wine to calibrate the disk or the calculation software and initialize the atmospheric pressure sensor. This operation must be repeated each time the device is turned on.

### • Rinsing

- Check that the drain valve is correctly closed.
- Refill the boiler with the standard wine to the bottom of the condenser. Then empty the boiler using the drain valve.

Never turn on the instrument using the heating switch when the boiler is empty or when draining. You might damage the instrument if you do not observe this rule.

### • Boiling point of the standard wine

- Refill the boiler to the bottom of the meniscus above the red line.
- Turn on the instrument using the on/off switch.

The main screen appears.

- Let cold water circulate into the condenser, with a regular flow, not excessive.
- Press the heating button during at least 8 seconds to launch calibration.

The following screen appears successively :

VERIFIER CHAUDIERE  
1111

CALIBRATION  
6666

For more accuracy, we advise you to turn on the instrument about 15 minutes before use.

The heating (shown by the symbol ) lasts about 7 minutes, during which the temperature rise can be seen on the screen.

To stop heating during the operation, press again on the heating button.

- When the boiling temperature is stable, the device makes a beep to indicate the end of the measurement. Heating stops automatically.
- The boiling temperature  $T_{\text{Standard}}$  remains on the display.

CALIBRATION  
93,10°C

The accuracy can be reached with standard wine, whose properties are similar to the wine to be tested. If no standard wine is available, it can be replaced by demineralized water déminéralisée, considered as standard wine with alcoholic content of 0% but accuracy will be lower.

- Match the alcohol content of the standard wine with its boiling temperature  $T_{\text{Standard}}$  on the disc, or use the programme «EBULLIOLOG», entering the values «Reference Temperature» and «Alcoholic reference».


- Drain.
- Press briefly the heating button to come back to the main screen. The device is ready for a sample measurement.

### • Example

Standard wine with official alcohol degree = 9,2 % vol ( $\pm 0,05\%$  vol.)

$T_{\text{Standard}} = 93,1^\circ\text{C}$

Calibrate the disc as indicated in the picture opposite.

*Model 160350T with traditional thermometer :*


*Press the heating button 5 seconds to launch heating. When the wine is boiling, check on the thermometer that the temperature is stable (after 1 to 2 minutes). Once it is stable, record the boiling temperature  $T_{\text{Standard}}$ . Turn off the heating using the heating switch. Adjust the disk or the programme as indicated above. Repeat calibration every 2 hours.*

## 4 - Sample measurement

A sample mesure can only be done after a calibration has been done, at least after turning the device on to reinitialize the atmospheric pressure sensor.

- **Rinsing**

- Check that the drain valve is correctly closed.
- Refill the boiler with the standard wine to the first ball of the condenser, then empty the boiler using the drain valve.

- **Boiling point of the standard wine**

- Refill the boiler to the bottom of the meniscus above the red line.
- Check that the refrigerating circuit is still active.
- Press the heating button during 5 to 8 seconds to launch measurement.

The heating (shown by the symbol ) lasts about 7 minutes, during which the temperature rise can be seen on the screen.

To stop heating during the operation, press again on the heating button.

- When the boiling temperature is stable, the device makes a beep to indicate the end of the measurement. Heating stops automatically.
- The boiling temperature  $T_{\text{Wine}}$  remains on the display.

- Check that the disc is correctly calibrated or the programme «EBULLIOLOG» correctly set up. Then read the alcohol content of the wine matching with the boiling temperature  $T_{\text{Wine}}$ , or in the programme «EBULLIOLOG», enter the value «Sample temperature» and click on the button Calculation of dosage to obtain the alcohol content of the sample.
- Drain.
- Press briefly the heating button to come back to the main screen. The device is ready for a new sample measurement.

- **Exemple**

Calibration of the disc with a standard wine at 9,2 % vol. and  $T_{\text{Standard}}=93,1^{\circ}\text{C}$  according to the previous example.

$T_{\text{Wine}} = 91,4^{\circ}\text{C}$

Alcohol content of the wine = 12,2% vol ( $\pm 0,1\%$ )

**Automatic pressure compensation**


After the calibration, an integrated atmospheric pressure sensor determines the correction to apply to the measured temperature. The temperature displayed at the end of the measurement is corrected yet from the possible atmospheric pressure changes bewteen calibration and measurement.

**Some wines with a high volume of CO<sub>2</sub> foam during boiling.**

Excessive foaming may interfere on the measure. Pour two drops of defoamer if necessary with the supplied dropping-pipette directly into the condenser to lessen foaming.

**Model 160350T with traditional thermometer :**

When the wine is boiling, check on the thermometer that the temperature is stable (after 1 to 2 minutes). Once it is stable, record the boiling temperature  $T_{\text{Wine}}$  and turn off the heating using the heating switch. Adjust the disk or the programme as indicated above.

## 5 - Maintenance

The boiler, heating cartridge and probe or thermometer must be cleaned once a month, or more frequently in case of intensive use.

- Drain the boiler and then close the drain valve.
- Leave the boiler and the probe (or thermometer) on the stand.
- Refill with cold water till the first ball in the condenser.
- Drain and then refill with the cleaning solution (Ref. 972500) to the angle of the boiler.
- Wait for 3 hours, or more in case of frequent use.
- Drain and then rinse thoroughly four or five times with water, until the product has completely disappeared.

THIS SOLUTION MUST BE USED BETWEEN 20 AND 50°C.

DO NOT TURN ON THE HEATING DURING THE CLEANING PART.

The Material Safety Data Sheet is available on request.

- *Never mouth-pipette this glass cleaner.*
- *The glass cleaner can cause burns. Rinse immediately down if the cleaning solution comes in contact with skin or eyes.*
- *The solution corrodes nonferrous metals. Do not use it to clean other equipments. Wear protective gloves and eye protection.*
- *The cleaning solution is a single-use product. Dilute and neutralize with a weak acid in case of throwing down the drain.*

## 6 - Warranty and return

- **Instruments under warranty should be returned:**

- In their original packaging (parcels and styrofoam);
- With shipping costs paid, including insurance against any breakage during the carriage.  
We do not warranty any failure due to an ordinary wear, a non-approved use, an on-site service on the instruments done by a third After Sales Service or to an external accident.

- **Instruments out of warranty should be returned:**

- In their original packaging when possible or perfectly secured in a double packing to avoid any direct contact between the different elements;
- With shipping costs paid, including insurance against any breakage during the carriage.

- **The warranty applies to :**

- All parts except the glassware, rubber hoses and the electronic probe, within the limits of a normal use of the instrument, as described in this user manual.

## 7 - Technical characteristics

| Ebulliometer with traditional thermometer - Ref. 160350T<br>Ebulliometer with electronic probe - Ref. 160350D |                                               |
|---------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| Sample type                                                                                                   | Dry wines                                     |
| Measure scale                                                                                                 | From 0 to 17 % Vol. Alcohol                   |
| Accuracy                                                                                                      | 0,1 % Vol. Alcohol                            |
| Dimensions and weight                                                                                         | 42 x 24 x 22 cm ; 3,2 kg                      |
| Input                                                                                                         | External supply 100-240 V - 47/63 Hz - 0,7 kg |
| Output                                                                                                        | 24 V DC                                       |
| Mercury thermometer                                                                                           | Scale 84 - 101°C / Accuracy 0,1°C             |
| Thermometer with probe                                                                                        | Type PT100 / Accuracy 0,1°C                   |

## 8 - Statement of conformity

**LABORATOIRES DUJARDIN-SALLERON**  
**872, route de la Gare**  
**37210 NOIZAY**  
**FRANCE**  
Tel : 00 33 2 47 25 58 25  
Fax : 00 33 2 47 25 58 30  
E-mail : info@dujardin-salleron.com

**STATEMENT OF CONFORMITY** 

*Type of instrument: Electric ebulliometer*

Model EBULLIOMETER Ref.160350T-160350D

Has been developed, designed and produced according to the regulations of the EMC Directive 2004/108/CE and of the Low Voltage Directive 2006/95/CE.

In accordance with the harmonized standards:

|                                   |                     |
|-----------------------------------|---------------------|
| EN 55022:2006                     | IEC 61000-4-4:2004  |
| EN 61000-3-2:2006                 | IEC 61000-4-11:2004 |
| EN 61000-3-3:1995+A1:2001+A2:2005 | IEC 61000-4-3:2006  |
| EN 55024:1998+A1:2001+A2:2003     | IEC 61000-4-5:2005  |
| EN 61010-1:2004                   |                     |

Arcueil, December 2009

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# 9 - Appendix

- **Optional accessories**

- |          |                                                    |
|----------|----------------------------------------------------|
| 160351   | Method for mellow wines                            |
| 160352   | Method for mellow wines and musts                  |
| 160353   | Method for ciders and pommeau                      |
| 160354   | Method for beers with low fermentation             |
| 160355/1 | Method for vinegar                                 |
| 160356   | Clé USB 1 Go EBULLIOLOG with calculation programme |

- **User manual for EBULLIOLOG (optional)**


- Doble-click on the icon of the file EBULLIOLOG to open the programme.

Set up :

- Click on the flag to choose your language.
- Click on the thermometer to obtain the right temperature scale.

Use :

- After calibration, enter the following values :
  - **Reference Temperature** : It is the boiling temperature of the standard wine.
  - **Alcoholic reference** : It is the alcoholic content of the standard wine.
- After measuring the boiling point of the sample, enter the value **Sample temperature**.
- Click on Calculation of dosage to obtain the **Alcohol content of the sample**.

- **Process for standard wine preparation**

- Prepare 10 liters of wine, stabilized and homogenized.
- Pour it in small flasks of 250 ml.
- Have one or two of these flasks checked at one/two laboratories with official method (distillation + densimetry) to obtain the value of the alcoholic degree.
- Keep the flasks protected from light, between 5°C and 25°C.

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